



Energy Transfer Solar

Model No. ET-8593




Table of Contents

Equipment List	3
Introduction	4
Equipment Description	4
Equipment Setup	5-6
Suggested Experiments	7
Experiment 1: Solar Heating and the Greenhouse Effect	7-8
Experiment 2: Solar Constant (Advanced Lab)	9-10
Sample Data/Results	11
Appendix A: Temperature/Resistance Conversion Table	12
Appendix B: Technical Support	13
Appendix C: Copyright and Warranty Information	13

Energy Transfer - Solar

Model No. ET-8593

Equipment List



Included Equipment	Replacement Model Number*
1. Solar Box , 10.75 x 8.25 x 2.50"	648-08412
2. Plastic Cover, 11.0 x 8.25 x 1.25"	650-065
3. Aluminum Plate, 6.50 x 9.0", 85 g	648-08413
4. Cable assembly	514-08366

*Use Replacement Model Numbers to expedite replacement orders.

Additional Equipment Required	
PASPORT™ Xplorer or a laptop computer	PS-2000
DataStudio® software	Various (See PASCO catalog)
Temperature Sensor or Thermistor Sensor or Ohmmeter/Multimeter	PS-2125 or CI-6527A (For ohmmeter, see PASCO catalog.)
Piece of cardboard (1 ft. square)	NA
A computer	NA

Introduction

The Energy Transfer-Solar box (ET-8593) can be used for demonstrating the concept of solar heating, including the greenhouse effect.

Equipment Description

a) Plastic Cover

The clear, plastic cover snaps onto the Solar box and acts as insulator to isolate and trap air inside, reduce convection currents, and demonstrate the greenhouse effect. The cover is very transparent to visible light but not infrared light.

b) Aluminum Plate

The aluminum plate is painted a non-reflective flat black that absorbs light very well. The hot aluminum plate re-radiates in the far infrared region, and thus the heat energy is trapped under the cover.

The reverse side of the aluminum plate is not painted. The plate can be flipped inside the box to study differences in solar heating and/or cooling between the aluminum and black surfaces. The aluminum plate can be removed to measure its mass. The white, plastic knob also serves as an indicator for the sun's angle. When the sun is perpendicular to the aluminum plate, no indicator shadow appears on the plate.

c) Solar Box

The Solar Box holds the aluminum plate and plastic cover. On the side of the box is a rod clamp for mounting the box to a rod stand. When mounted to a rod stand, the box can be adjusted to the sun's angle.

d) Thermistor

Inside the Solar Box is a 10K thermistor for measuring temperature. The thermistor cables are not removable from the box. The thermistor contact (metal lug) is fastened in the center, on the underside of the aluminum plate. The side jacks on the Solar Box allow you to connect a Temperature Sensor or ohmmeter to the thermistor.



Note: The thermistor's temperature range is -35°C to $+135^{\circ}\text{C}$ (242,800 to 265 ohms).

Equipment Setup

Note: A rod stand is required for adjusting the Solar Box to the sun's angle.

1. Remove the thumbscrew on the aluminum plate. Place the thermistor lug underneath the center hole on the aluminum side of the plate. Insert the thumbscrew through the hole. On the black side, put the shadow indicator over the screw and tighten (See Figure 1).




Figure 1: Thermistor position

2. Place the aluminum plate inside the Solar Box, with the black side face up. (**Note:** Keep the thermistor contact on the underside of the plate.)
3. Snap the bottom and top tabs of the clear, plastic cover onto the Solar Box. (See Figure 3).

4. Using the thermistor cable, connect the Temperature Sensor (or Thermistor Sensor) to the two jacks on the side of the Solar Box (Figure 2). If a Temperature Sensor is not available, connect an ohmmeter to the side jacks.




Figure 2: Connecting the thermistor cable

5. Use the rod clamp (on the side of the Solar Box) to mount the base of the Solar Box to a rod stand (Figure 3).
6. Adjust the angle of the box such that the sun's rays enter the box perpendicularly. Use the white knob indicator as a guide.



Figure 3: Mounting the Solar Box to a rod stand

Note: If there is no shadow on the plate, the sun's rays are perpendicular to the plate.

7. Plug the Temperature Sensor into a PASPORT Xplorer. To take a temperature reading, click the **Start** button on Xplorer. (**Note:** The Xplorer data can later be uploaded into DataStudio and viewed in a DataStudio graph display.)

OR

If you are using an ohmmeter (instead of a Temperature Sensor), turn on the meter and take a resistance measurement. To find the temperature, use the resistance-to-temperature conversion chart in Appendix A.).

OR

If you are using a laptop and a temperature probe, plug the Temperature Sensor plug into a USB port on your laptop computer. Launch DataStudio and click the **Start** button to collect data.



WARNING: To avoid burns or bodily injury, when heating the box, do not overheat the box (above 100°C) and do not touch either side of the aluminum plate or the thermistor contact.



CAUTION: Overheating the box may permanently damage the thermistor and the plastic lid. The thermistor's maximum temperature capacity is 135°C.

Suggested Experiments

Experiment 1: Solar Heating and the Greenhouse Effect

Equipment Required	
Energy Transfer - Solar (ET-8593)	Temperature Sensor (PS-2125) or Thermistor Sensor (CI-6527A) or an Ohmmeter/Multimeter
Rod Stand (ME-9355)	PASPORT Xplorer (PS-2000) or laptop computer
Piece of cardboard	Temperature vs. Resistance Chart (See Appendix A)
DataStudio Software	

Part I: Solar Heating

1. Mount the box with plate to a rod stand, such that the Sun's angle is perpendicular to the aluminum plate and the white plastic knob has no shadow. Keep the black side of the aluminum plate facing up (See Figure 4).
2. Use a PS-2000 Xplorer or a laptop computer for data collection. Set the sample rate in either Xplorer or DataStudio for 2 Hz.
3. Have a piece of cardboard available to shade the box while setting up.
4. In DataStudio, click the **Start** button to begin data collection and remove the cardboard shade.
5. With the plastic cover on, take a run of data in DataStudio. Let the box heat until the temperature levels off. (The approximate duration is 10 to 30 minutes, depending on the outside temperature and the intensity of the sunlight.)



Figure 4: Setup for Solar Heating Experiment



CAUTION: Overheating the box may permanently damage the thermistor and the plastic lid. The thermistor's maximum temperature capacity is 135°C.

Note: Watch the angle of the sun. *The angle of the sun must be 90 degrees to the box while you are collecting data.* You might have to adjust the angle of the box during the run.

6. Repeat step 5 with the plastic cover off.

Analysis

1. Look carefully at both curves at the start of the run. The slope (rate of heating) for the uncovered box should be larger than for the covered box. Why?
2. Which has the highest final temperature, the covered box or the uncovered box?
3. Which curve has a more constant heating rate? Why?

Part II. Solar Heating Comparison: Aluminum vs. Black Surface

Compare the aluminum side up to black side up with the cover on. Which surface is a better absorber of energy? Look at not only how fast the plate heats up, but collect data long enough to look at the final temperature. The black side should heat up much faster than the aluminum side, but does the black side reach a higher final temperature?

Part III. Cooling Comparison: Aluminum vs. Black Surface

For both sides (aluminum and black), start with the plate hot (Let it sit in the sun), and then move the plate to the shade to watch it cool. Which surface cools faster? Which is a better emitter of energy? Try cooling both with and without the cover on the Solar Box.



WARNING: To avoid burns or bodily injury, when heating the box, do not overheat the box (above 100°C) and do not touch either side of the aluminum plate or the thermistor contact.

Experiment 2: Solar Constant

Equipment Required	
Energy Transfer - Solar (ET-8593)	Temperature Sensor (PS-2125) or Thermistor Sensor (CI-6527A)
Rod Stand (ME-9355)	PASPORT Xplorer (PS-2000) or laptop computer
Piece of cardboard	Measuring tape and scale
DataStudio Software	

***Note:** This is a more advanced lab. Two Temperature Sensors or one temperature sensor and a thermometer are required. You will use one Temperature Sensor to measure the temperature of the aluminum plate and a second Temperature Sensor (or thermometer) to measure the ambient temperature.

1. Disconnect the thermistor from the plate and measure the mass of the plate. Measure the plate's size and calculate the area of plate.
2. Cool the plate to 10°C to 20°C below the outside temperature. (You can stick the plate in a refrigerator or use an ice cube). Be sure the plate is dry.
3. Place the aluminum plate in the Solar Box with the black side facing up to the sun. **Do not use the plastic cover.**
4. Connect the Temperature Sensor to the side jacks of the box with the supplied cable. (If possible, have a second temperature sensor measuring outside temperature. **Note:** The second Temperature Sensor (or a thermometer) must be in the shade for an accurate determination of the outside ambient air temperature.)
5. Recheck the sun's angle. You might have to adjust the box relative to the sun's angle during the run. (**Note:** The angle of the sun relative to the box must be 90 degrees.)
6. In DataStudio, create a graph of temperature vs. time. For the time variable, use seconds (not minutes) on the graph.
7. In DataStudio, click the **Start** button to begin recording (at the default sample rate of 2 Hz.). Heat the box until it is 10°C to 20°C above outside temperature.



WARNING: To avoid burns or bodily injury, when heating the box, do not overheat the box (above 100°C) and do not touch either side of the aluminum plate or the thermistor contact.



CAUTION: Overheating the box may permanently damage the thermistor and the plastic lid. The thermistor's maximum temperature capacity is 135°C.

Analysis

a) Find the slope of the line tangent to the curve at outside ambient air temperature. Do not use the Slope Tool. Highlight a small section at outside ambient temperature. At ambient temperature, the heating is only being caused by the sunlight. Below the ambient temperature, the surrounding air is *cooling* the container. Above the ambient temperature, the surrounding air is *heating* the container.

$$\text{slope} = \frac{\Delta T}{\Delta t} = \text{change in temperature/change in time}$$

b) Theory for heat flow

$Q = mc\Delta T$ where Q =heat, c =specific heat, and ΔT = change in temperature

$$\text{Power} = \frac{Q}{\Delta t} = mc \frac{\Delta T}{\Delta t}$$

$$\text{Intensity} = \frac{\text{Power}}{\text{Area}} = \frac{mc \frac{\Delta T}{\Delta t}}{\text{Area}} \quad \text{where } \frac{\Delta T}{\Delta t} \text{ is the slope of the graph.}$$

c) Using your slope and the other measured quantities, calculate the intensity of the sun's light. The intensity (solar constant) at the top of the Earth's atmosphere is about 1400 Watts/m². On a good clear day with the sun high in the sky, you can get over 1000 Watts/m² on the surface.


For Further Study:

- Compare the intensity at noon to later in the day.
- Compare a clear day to a slightly overcast day.
- Compare summer to winter.


Sample Data/Results

Experiment 1: Solar Heating

The covered box has the final higher temperature. With the covered box, the greenhouse effect occurs. Light enters the transparent cover, but the infrared light is not re-radiated back out. The infrared light heats the plate. Also, the cover traps the hot air inside the box, isolating and insulating the air.



Experiment 2: Solar Constant



Note: This data was taken late in the afternoon.

$$\text{Power} = \frac{Q}{\Delta t} = mc \frac{\Delta T}{\Delta t} = (0.085 \text{ kg}) (900 \text{ joules/kg}^\circ\text{C}) (0.361) = 27.6 \text{ watts}$$

$$\text{Intensity} = 27.6 \text{ watts}/0.038 \text{ m}^2 = 730 \text{ watts/m}^2$$

Appendix A: Resistance/Temperature Conversion Table

Resistance (Ohms)	Temperature (Celsius)	Resistance (Ohms)	Temperature (Celsius)	Resistance (Ohms)	Temperature (Celsius)
32,660	0	6,808	34	1,876	68
31,040	1	6,532	35	1,813	69
29,500	2	6,268	36	1,751	70
28,060	3	6,016	37	1,693	71
26,680	4	5,776	38	1,637	72
25,400	5	5,546	39	1,582	73
24,180	6	5,326	40	1,530	74
23,020	7	5,118	41	1,480	75
21,920	8	4,918	42	1,432	76
20,880	9	4,726	43	1,385	77
19,900	10	4,544	44	1,341	78
18,970	11	4,368	45	1,298	79
18,090	12	4,202	46	1,256	80
17,260	13	4,042	47	1,216	81
16,460	14	3,888	48	1,178	82
15,710	15	3,742	49	1,141	83
15,000	16	3,602	50	1,105	84
14,320	17	3,468	51	1,071	85
13,680	18	3,340	52	1,038	86
13,070	19	3,216	53	1,006	87
12,490	20	3,098	54	975	88
11,940	21	2,986	55	945	89
11,420	22	2,878	56	916	90
10,920	23	2,774	57	889	91
10,450	24	2,674	58	862	92
10,000	25	2,580	59	836	93
9,574	26	2,488	60	811	94
9,166	27	2,400	61	787	95
8,778	28	2,316	62	764	96
8,408	29	2,234	63	742	97
8,058	30	2,158	64	720	98
7,722	31	2,082	65	699	99
7,404	32	2,012	66	679	100
7,098	33	1,942	67		

Appendix B: Technical Support

For assistance with the ET-8593 Energy Transfer - Solar or any other PASCO products, contact PASCO as follows:

Address: PASCO scientific

10101 Foothills Blvd.

Roseville, CA 95747-7100

Phone: (916) 786-3800

FAX: (916) 786-3292

Web: www.pasco.com

Email: techsupp@pasco.com

Appendix C: Copyright and Warranty Information

Copyright Notice

The PASCO scientific 012-08428A *Energy Transfer - Solar Manual* is copyrighted and all rights reserved. However, permission is granted to non-profit educational institutions for reproduction of any part of the 012-08428A *Energy Transfer - Solar Manual* providing the reproductions are used only for their laboratories and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Limited Warranty

PASCO scientific warrants the product to be free from defects in materials and workmanship for a period of one year from the date of shipment to the customer. PASCO will repair or replace, at its option, any part of the product which is deemed to be defective in material or workmanship. The warranty does not cover damage to the product caused by abuse or improper use. Determination of whether a product failure is the result of a manufacturing defect or improper use by the customer shall be made solely by PASCO scientific. Responsibility for the return of equipment for warranty repair belongs to the customer. Equipment must be properly packed to prevent damage and shipped postage or freight prepaid. (Damage caused by improper packing of the equipment for return shipment will not be covered by the warranty.) Shipping costs for returning the equipment after repair will be paid by PASCO scientific.